



Horizon® - Superior image quality and value for CR/DR

Satisfy Radiologists

- **Print Versatility-** Horizon incorporates film and paper into one incredibly compact and versatile medical imager. Horizon offers several models specifically designed for the needs of CR/DR. Horizon GS prints 14" x 17" film and low cost DirectVista® Grayscale Paper while Horizon GS-s prints 11" x 14" film and grayscale paper. Both models also print 8" x 10" film, and A/A4 grayscale paper, allowing you to choose the media size that fits your image, reducing cost and waste.

Lower Total Cost

- **Calculate Savings-** Horizon provides lower hardware, film and service costs than any competing dry imager. Horizon also has the industry's best warranty and first year service is included.
- **Low Cost Media -** Codonics DirectVista Paper is a high quality media ideal for referring physicians at a fraction of the cost of film available in 14" x 17" and 11" x 14" sizes, providing a cost-saving alternative to large format film. Additionally A/A4 letter size paper offers a **savings of up to 70% over film.**

Improve Workflow

- **Speed -** With the fastest time-to-first-print, Horizon prints up to 100 films an hour.





Horizon®

Superior image quality and value for CR/DR.

Horizon instantly delivers diagnostic film and cost-saving paper, allows you to load three types of media at the same time and holds up to 300 sheets of film. This unprecedented variety in printing options creates a new freedom of choice and allows users to match the output with the need, resulting in solutions capable of immediately reducing costs, minimizing waste and maximizing workflow. Horizon packs this all into a compact device, weighing less than 70 pounds and taking up less than two feet of counter space, eliminating costly siting requirements.




Radiology / Imaging Center

- Direct or Network Connection
- Compact and Easy to Site
- Up to 24 Simultaneous DICOM Connections
- Satisfy Demands from Referring Physicians
- High Throughput
- Reduce Costs

Referring Physicians

- Clinical Review
- Patient Consults
- Patient Medical Files
- Surgical Planning
- Patient Marketing
- Ideal for Orthopaedic True-Size Images



For more information contact Codonics at **800-444-1198** or **440-243-1198** and ask for our CR/DR Global Sales Manager.

Copyright © 2005 Codonics, Inc.



17991 Englewood Drive
 Middleburg Heights, OH 44130 USA
 (440) 243-1198
 (440) 243-1334 Fax
 Email info@codonics.com
 www.codonics.com

900-244-001 Rev B