

Horizon® Imager Software

Release Notes

Version 1.7.1 Release Notes

With version 1.7.0 software, certain rare types of print engine errors required the system to be rebooted to restore normal operation. Version 1.7.1 software repairs this defect by more thoroughly re-initializing the print engine when these types of errors occur, such that a reboot is not required.

In all other respects, version 1.7.1 is identical to version 1.7.0. See Version 1.7.0 Horizon Imager Software Release Notes (Codonics Part Number 900-219-001) for full details.

This software is compatible with all Horizon imagers. Any Horizon imager running software versions 1.4.x through 1.7.x can be upgraded to version 1.7.1.

Get it all with just one call
1-800-444-1198

All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice. Patents pending.

Copyright © 2000-2003 by Codonics, Inc. Printed in the U.S.A. Part No. 901-084-002 Rev.A.



We bring the future into focus

17991 Englewood Drive
Middleburg Heights, OH 44130 USA
(440) 243-1198
(440) 243-1334 Fax
Email info@codonics.com
www.codonics.com

Codonics Limited KK
New Shibaura Bldg. F1
1-3-11, Shibaura
Minato-ku, Tokyo, 105-0023 JAPAN
Phone: 81-3-5730-2297
Fax: 81-3-5730-2295