

Horizon Imager Software

Release Notes

Version 1.8.7 Release Notes

(Changes since Version 1.8.6)

Horizon® 1.8.7 software allows Horizon to support hard drives larger than 34 GB.

This software is compatible with all Horizon imagers including LifeImager products. Any Horizon imager running software versions 1.5.x through 1.8.6 can be upgraded to version 1.8.7. Imagers running version 1.4.x software must first be upgraded to version 1.6.0 before loading version 1.8.7. If your Horizon imager is running version 1.4.x software, please contact Codonics Technical Support.

This document details only changes from versions 1.8.6 software. If you are upgrading from an older software version, previous Horizon Imager Software Release Notes can be downloaded from the Codonics web site.

Product Improvements

This section details features which have been improved with Horizon version 1.8.7 software.

Horizon 1.8.7 software patches the main operating system, Red Hat Linux 7.3, kernel version 2.2.5, to support hard drives larger than 34 GB.

Get it all with just one call
1-800-444-1198

All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice. Patents pending.

Copyright © 2008 by Codonics, Inc. Printed in the U.S.A. Part No. 901-084-013 Rev. 01.



17991 Englewood Drive
Middleburg Heights, OH 44130 USA
(440) 243-1198
(440) 243-1334 Fax
Email info@codonics.com
www.codonics.com

Codonics Limited KK
New Shibaura Bldg. F1
1-3-11, Shibaura
Minato-ku, Tokyo, 105-0023 JAPAN
Phone: 81-3-5730-2297
Fax: 81-3-5730-2295