

# Codonics® Virtua™ Operating Software

**Release Notes**

---

## **Version 1.2.4 Release Notes**

### **Summary**

Version 1.2.4 software adds support for a Japanese User Interface. This software is compatible with all Virtua systems.

Loading 1.2.4 software will change the User Interface from English to Japanese provided that the locale is set for Japan (e.g. ja\_JP). If a customer wishes to return to English or to change the locale after loading 1.2.4 software, call Codonics technical support.

### **New Features**

- ◆ **Support for Japanese Language User Interface.** The Virtua user interface including labels, button text, dialog messages, column headings, and warning messages are translated into Japanese. The translation is a combination of Kanji, Hiragana, and Katakana. There are a few exceptions (job profile names, label templates, and disc features) which remain in English.

### **Technical Support**

If problems occur during software installation, contact Codonics Technical Support between the hours of 8:30AM and 5:30PM EST (Weekends and U.S. Holidays excluded).

Phone: 440-243-1198  
Email: support@codonics.com  
Website: www.codonics.com

---

*Get it all with just one call*  
**1-800-444-1198**

All registered and unregistered trademarks are the property of their respective owners. Specifications subject to change without notice.

Copyright © 2005,2007 by Codonics, Inc. Printed in the U.S.A. Part No. 901-150-006 Rev. A.



*We bring the future into focus*

17991 Englewood Drive  
Middleburg Heights, OH 44130 USA  
(440) 243-1198  
(440) 243-1334 Fax  
Email info@codonics.com  
www.codonics.com

Codonics Limited KK  
New Shibaura Bldg. F1  
1-3-11, Shibaura  
Minato-ku, Tokyo, 105-0023 JAPAN  
Phone: 81-3-5730-2297  
Fax: 81-3-5730-2295